


# Here is the CENTRISCAN 75

A simpler and less expensive version of the classical analytical ultracentrifuge — but with a built-in scanner and full preparative facilities


#### Analytical

*Absorption scanning down to 254 nm  
Refractometric scanning technique  
giving Schlieren traces*



#### Preparative

*Preparative rotors  
to 75,000 r.p.m.  
500,000 x 'g'*



The MSE Centriscan 75 is no compromise; it combines full preparative and analytical facilities — one centrifuge to buy instead of two. On the analytical side it is capable of work of the highest precision and sensitivity, recording absorbance of up to four samples in the same rotor.

The optical system is a unique form of split-beam scanner operating down to 254nm and scans are drawn directly on a recorder — no cameras! A simple changeover converts the instrument to refractometric scanning and permits the recording of Schlieren diagrams. On the preparative side, 75,000 rpm gives

500,000 x 'g' — and unique versatility. Ask now for Publication No. 501

**Measuring & Scientific Equipment Ltd., Manor Royal, Crawley, Sussex.**


**Telephone: Crawley (STD 0293) 31100**





**GALACTOSE**

**GLUCURONIDES**



For the first time, all enzymes and substrates for the "Leloir Pathway", Galactose Metabolism and Glucuronide Synthesis are available. Many of the following preparations are offered here for the first time anywhere. Prices are tentative. Many have already been substantially reduced. Others will be reduced as demand warrants. Since most enzymes are about 30% more active at 30°C than at 25°C, reactions performed at higher temperatures may result in saving time and material.

**GALACTOKINASE**

G0130 From Yeast. Present lot contains 75% Buffer Salts.

Approx. 30-50 units/mg Protein

Unit Definition: One unit will convert one  $\mu$ Mole of Galactose to Gal-1-P

per min. at pH 7.0 at 25°C.

Approx. 1% Hexokinase impurity

2 units	\$ 5.50	20 units	\$25.00
5 units	10.00	100 units	75.00

**GALACTOSE-1-PHOSPHATE**

**URIDYL TRANSFERASE**

G8128 From Yeast. Present lot contains approx. 85% Buffer

Salts. Approx. 2 units/mg Protein

Unit Definition: One unit will form one  $\mu$ Mole of G-1-P from UDPG and

Gal-1-P per min. at pH 8.7 at 25°C.

Approx. 5% UDPG Pyrophosphorylase impurity

2 units	\$5.00	50 units	\$75.00
---------	--------	----------	---------

**PHOSPHOGLUCOMUTASE**

P7502 From Rabbit Muscle. Crystalline. Ammonium Sulfate

Suspension. 80-130 units/mg Protein.

Unit Definition: One unit will convert one  $\mu$ Mole G-1-P to G-6-P per

min. at 7.4 at 30°C.

200 units	\$ 7.90	1000 units	\$28.00
500 units	18.00	2500 units	56.00

**URIDINE-5'-DIPHOSPHOGALACTOSE-**

**4-EPIMERASE**

U3251 From Yeast. Present lot contains approx. 40% Buffer

Salts. 10-15 units/mg Protein

Unit Definition: One unit will convert one  $\mu$ Mole UDPGal to UDPG per

min. at 25°C, pH 8.8.

Approx. 0.5% Galactokinase impurity

5 units	\$5.00	20 units	\$12.00
		100 units	\$50.00

Also available—See our Catalog or inquire:

- D(+)-GALACTOSE
- $\alpha$ -D-GALACTOSE-1-PHOSPHATE (synthetic)
- $\alpha$ -D-GLUCOSE-1-PHOSPHATE
- D-GLUCOSE-6-PHOSPHATE

**URIDINE-5'-DIPHOSPHOGLUCOSE**

**DEHYDROGENASE**

U5500 Powder. Buffered preparation. From Bovine liver.

Approx. 500 units/mg Protein

Unit Definition: One unit will oxidize  $4.0 \times 10^{-5}$   $\mu$ Moles of UDPG/min.

at pH 8.7 at 25°C.

10,000 units	\$9.90	20,000 units	\$16.50
		50,000 units	\$33.00

**URIDINE-5'-DIPHOSPHOGLUCOSE**

**PYROPHOSPHORYLASE**

U1501 From Yeast. Ammonium Sulfate Suspension. Ap-

prox. 45 units/mg Protein

Unit Definition: One unit will cause the formation of 1.0  $\mu$ Mole of

G-1-P from UDPG and Inorganic Pyrophosphate per min. at pH 7.6 at

25°C.

Impurities include Alcohol Dehydrogenase and G6PD

25 units	\$ 6.50	500 units	\$ 60.00
100 units	18.00	1000 units	100.00

**URIDINE-5'-DIPHOSPHOGLUCURONYL**

**TRANSFERASE**

U3626 From Rabbit liver. Crude microsomal preparation.

Present lot contains approx. 10% Buffer Salts. Approx. .0015

units/mg Protein

Unit Definition: One unit will transfer 1.0  $\mu$ Mole of Glucuronic Acid

from UDPGA to Phenolphthalein per min. at pH 8.0 at 37°C.

Present lots contain Glucuronidase impurity approx. 0.001  $\mu$ Molar units

per mg at pH 4.5; approx. 0.00005  $\mu$ Molar units/mg at pH 8.0.

5 units	\$6.50	25 units	\$30.00
---------	--------	----------	---------

**$\alpha$ -D-GALACTOSE-1-PHOSPHATE**

G0380 Dipotassium Salt. Enzymatically prepared. Purity:

98-100%. Substantially free of G-1-P.

10 mg	\$ 4.80	250 mg	\$ 54.00
50 mg	16.20	500 mg	90.00
100 mg	27.00	1 g	150.00

100 mg	27.00	1 g	150.00
--------	-------	-----	--------

SEVERAL GLUCURONIDES

- URIDINE-5'-DIPHOSPHOGALACTOSE
- URIDINE-5'-DIPHOSPHOGLUCOSE
- URIDINE-5'-DIPHOSPHOGLUCURONIC ACID
- URIDINE-5'-TRIPHOSPHATE

**ORDER DIRECT - TELEPHONE COLLECT from ANYWHERE in the WORLD**  
 Day, Station to Station, 314-771-5750 — Night, Person to Person, Dan Broida, 314-993-6418  
 TWX (Teletype) Day or Night: COLLECT 910-761-0593  
 TELEGRAM: SIGMACHEM, St. Louis, Missouri


The Research Laboratories of  
**SIGMA** CHEMICAL COMPANY

MAILING ADDRESS: P. O. BOX 14508, ST. LOUIS, MO., 63178, U.S.A.

MANUFACTURERS OF THE FINEST BIOCHEMICALS AVAILABLE

Distributed through:  
 SIGMA LONDON Chem. Co. Ltd. • 12, Lattice St., London, S.W.6., England • Telephone: 01-736-5823 (Reverse Charges)  
 SIGMA ISRAEL Chem. Co. Ltd. • 28 Kat-Gimel St., Givataim, Israel • Telephone: 03 766654 (Reverse Charges)

# SYSTEMATICS



## A Systematic approach to chromatography All components for complete systems

Pharmacia Fine Chemicals offers a wide range of compatible components for a variety of gel filtration and ion exchange chromatographic systems.

- A** COLUMNS. A comprehensive range of jacketed and unjacketed columns available in diameters from 0.9 to 21.5 cm. Lengths from 15 to 100 cm. Special solvent resistant columns also available in diameter 2.5 cm. Lengths 45 and 100 cm.
- B** FLOW ADAPTORS. Permit upward flow elution, sample application, recycling, flow reversal, etc.
- C** VALVES. 3- and 4-way valves allow easy operation and flexibility for sample application, recycling, flow reversal, by-pass, etc.
- D** RESERVOIRS. Bi-functional gel and eluant reservoirs for column packing and elution at constant pressure.
- E** CONNECTORS. 2- and 3-port connectors for capillary tubing allow easy interconnection of chromatographic components without flanging or flaring.
- F** FLOW CONTROL VALVE. Designed for shut-off as well as flow control.

### NEW

Jacketed columns with new design. Diameters 1.6 and 2.6 cm (K16 and K26). Lengths 40, 70 and 100 cm. Special 1.6 x 20 cm column for ion exchange chromatography.

Flow adaptors for 1.6 and 2.6 cm columns (A 16 and A 26).

Flow control valve (FCV-1). Standard on 1.6 and 2.6 cm columns.

Pharmacia (Great Britain) Ltd  
Paramount House  
75 Uxbridge Road Ealing W5 5SS  
Telephone: 01-579 0102/7



**Pharmacia**  
**Fine Chemicals**

# NEW OXOID electrophoresis strips are now so good, we just had to put a brand- new container round them.

## New purer cellulose acetate permits—

- \* more clearly defined separations
- \* improved and permanent clearance of strips



## New Oxoid electrophoresis strips are now—

- \* interleaved for easier handling in a more practical container
- \* available in a wider range of sizes

For further information, write to  
Oxoid Limited, Southwark Bridge Road,  
London SE1 9HF. Telephone 01-928 4515.

**OXOID** for progressive laboratories

# Hoppe-Seyler's Zeitschrift für Physiologische Chemie

Editors in Chief

A. BUTENANDT · F. LYNEN · G. WEITZEL

Subscription Rates

For one volume (12 parts) DM 420,—

Vol. 352 No. 11

Contents

November 1971

- A DNA polymerase from *Ehrlich* ascites tumor cells with preference for native DNA**  
K. HARM, H. HILZ and M.-M. OLDEKOP
- Isolation and purification of a calcium-binding protein from human tissues**  
P. PIAZOLO, M. SCHLEYER and H. E. FRANZ
- Bis(t-butylloxycarbonyl)insulin**  
R. GEIGER, H.-H. SCHÖNE and W. PFAFF
- Studies in vitro on the biosynthesis of corticosteroids in rat adrenal gland and on the action of ACTH, II: Comparison of the synthesis from endogenous precursors with that from exogenous [4-<sup>14</sup>C]cholesterol, [4-<sup>14</sup>C]pregnenolone and [4-<sup>14</sup>C]progesterone**  
D. LOMMER, R. I. DORFMAN and E. FORCHIELLI
- The purification of thymidine kinase from *Escherichia coli* by affinity chromatography**  
W. ROHDE and A. G. LEZIU
- Inhibition of acylneuraminat-pyruvatylase: Evidence of intermediary Schiff base formation and of a possible role of histidine residues**  
R. SCHAUER and M. WEMBER
- The properties of carboxypeptidase C from orange leaves**  
B. SPRÖSSLER, H.-D. HEILMANN, E. GRAMPP and H. UHLIG
- Metabolism of sphingosine bases, XVII: Stereospecificities in the introduction of the 4t-double bond into sphinganine yielding 4t-sphingenine (sphingosine)**  
W. STOFFEL, G. ASSMANN and K. BISTER
- Ribosomal proteins, XI: Isolation of proteins from 70S ribosomes of *Escherichia coli***  
E. KALTSCHMIDT, V. RUDLOFF, H.-G. JANDA, M. ČECH, K. NIERHAUS and H. G. WITTMANN
- Chromophoric substrates, IV: A new substrate for the quantitative determination of chymotrypsin A**  
E. WÜNSCH, A. HÖGEL-BETZ and A. JAEGER
- Chromophoric substrates, V: The determination of collagenase with N-[(azulen-4-yl)acetyl]peptides**  
E. WÜNSCH, E. JAEGER and G. SCHÖNSTEINER-ALTMANN
- Chromophoric substrates, VI: On the specificity of collagenase**  
E. WÜNSCH, E. JAEGER and G. SCHÖNSTEINER-ALTMANN
- Chromophoric substrates, VII: On the specificity of carboxypeptidase B**  
E. WÜNSCH, E. JAEGER and G. SCHÖNSTEINER-ALTMANN
- Chromophoric substrates, VIII: Mass spectrometric investigations on N-[(azulen-4-yl)acetyl]peptides**  
E. WÜNSCH and E. JAEGER
- Protease inhibitors in human sperm plasma. Isolation by affinity chromatography and inhibition characteristics**  
E. FINK, E. JAUMANN, H. FRITZ, H. INGRISCH and E. WERLE
- [B1-p-iodophenylalanine]insulin, a homogeneous moniodoinsulin**  
G. KRAIL, D. BRANDENBURG, H. ZAHN and R. GEIGER
- Pteridines in the honey bee *Apis mellifica*: Isolation of new lumazines**  
J. H. DUSTMANN
- Chromatographic separation of gangliosides on precoated silicagel thin-layer plates**  
D. H. VAN DEN EIJNDEN



Walter de Gruyter · Berlin · New York

# CLINICAL SCIENCE

EDITORIAL BOARD

## *For the Biochemical Society*

R. HOFFENBERG, *Chairman*

C. N. HALES, R. G. HUNTSMAN, G. H. LATHE, J. LIDDELL,  
I. MACINTYRE, K. L. MANCHESTER, J. A. OWEN, R. W. E. WATTS

## *For the Medical Research Society*

J. S. ROBSON, *Deputy Chairman*

R. D. COHEN, W. I. CRANSTON, G. CUMMING, D. M. MATTHEWS,  
A. POLAK, J. I. S. ROBERTSON, B. ROBINSON, S. J. G. SEMPLE

VOLUME 42, No. 2

FEBRUARY 1972

### CONTENTS

Breathing patterns and regional ventilation distribution in tetraplegic patients and in normal subjects.  
By B. BAKE, A. R. FUGL-MFYER and G. GRIMBY

The effects of synthetic human and salmon calcitonins on electrolyte excretion in the rat. By C. C. WILLIAMS, E. W. MATTHEWS, JANE M. MOSELEY and I. MACINTYRE

Urinary excretion of nitrogen from [<sup>15</sup>N]valine, [<sup>15</sup>N]leucine and [<sup>15</sup>N]isoleucine in malnourished and recovered children. By W. W. C. READ, D. S. McLAREN and MARIE TCHALIAN

The antidiuretic action of diazoxide. By J. E. F. POHL, H. THURSTON and J. D. SWALES

Studies in man with <sup>18</sup>F. By D. J. HOSKING and M. J. CHAMBERLAIN

Search for evidence of lung chemoreflexes in man: study of respiratory and circulatory effects of phenyldiguamide and lobeline. By S. K. JAIN, S. SUBRAMANIAN, D. B. JULKA and A. GUZ

Electrolyte changes in patients with subarachnoid haemorrhage. By I. M. JAMES

The influence of numerous pregnancies and lactations on bone dimensions in South Africa Bantu and Caucasian mothers. By A. R. P. WALKER, BARBARA RICHARDSON and FAITH WALKER

The influence of disodium ethane-1-hydroxy-1,1-diphosphonate (EHDP) on the development of experimentally induced urinary stones in rats. By D. FRASER, R. G. G. RUSSELL, ORTRUN POHLER, W. G. ROBERTSON and H. FLEISCH

Orthostatic changes of haemodynamics, renal function, plasma catecholamines and plasma renin concentration in normal and hypertensive man. By M. MOLZAHN, TH. DISSMANN, S. HALIM, F. W. LOHMANN and W. OELKERS

Differential effect of noradrenaline and renal nerve stimulation on vascular resistance in the dog kidney and the release of a prostaglandin E-like substance. By J. C. MCGIFF, K. CROWSHAW, N. A. TERRAGNO, K. U. MALIK and A. J. LONIGRO

#### SHORT COMMUNICATION

A calcitonin-resistant primary action of parathyroid hormone in rats: calcium shift into bone. By C. J. ROBINSON, B. RAFFERTY and J. A. PARSONS

MEDICAL RESEARCH SOCIETY (Annual General Meeting held on 10 December 1971)

Demonstrations  
Symposium on the Red Cell  
Communications

*Subscription: £2.00 (\$7.50) per part; £20.00 (\$70.00) per year.*

*Orders may be placed with your bookseller, or sent direct to the publishers.*

**BLACKWELL SCIENTIFIC PUBLICATIONS LTD**

5 ALFRED STREET, OXFORD, OX1 4HB, ENGLAND

# Again our busy research workers bring you more new biochemicals



## 4 New Substrates

plus 3 other new  
biochemicals

## 9 New Molecular Biologicals

plus BD Cellulose.

## 20 New enzymes

plus two enzyme  
inhibitors.

## Combisets

Biochemical sets  
for standardisation  
and comparison.

### Enzymes

$\beta$ -amylase  
L-arginase  
Arylsulphatase  
Bromelain  
Carbonic Anhydrase  
Carboxypeptidase-A  
Carboxypeptidase-B  
Deoxyribonuclease Grade I  
Deoxyribonuclease Grade II  
Enzygel® - Papain  
 $\beta$ -glucosidase  
 $\beta$ -glucuronidase  
Neuraminidase  
Nuclease  
Nucleoside-diphosphate kinase  
Papain  
Pepsin  
Pepsinogen  
Phosphatase alkaline  
UDPG pyrophosphorylase

### Enzyme Inhibitors

Trypsin inhibitor-soya bean  
Trypsin inhibitor-egg white

### New Substrates

Ethanolamine  
 $\beta$ -glucose-6-phosphate  
DL-3-Hydroxybutyrate  
Succinyl-L-phenylalanine p nitranilide

### Molecular Biologicals

BD Cellulose  
Ribonucleic acid (soluble)

### Nucleosides

2'-deoxy-adenosine  
2'-deoxy-cytidine  
2'-deoxy-guanosine  
2'-deoxy-thymidine

### Nucleotides

Denry thymidine-5'-diphosphate  
Uridine-5'-diphospho-galactose  
Poly A-C  
Poly I-C  
Xanthosine 3':5'-monophosphate cyclic

### Combisets

Casein by Hammarsten  
Calibration Protein set  
Haemoglobin by Anson  
Insulin Chain A  
Insulin Chain B  
Insulin Chains A & B

### Other Biochemicals

Coenzyme A Grade II  
Insulin Chain A  
Insulin Chain B



For further details  
on any of these  
new products  
please contact:

Boehringer Corporation (London) Ltd.,  
Bilton House, Uxbridge Road,  
London, W.5.  
If you live in Britain or the Common-  
wealth

or  
Boehringer Mannheim G.m.b.H.,  
68, Mannheim 31,  
P.O.B. 51, West Germany  
if you live outside the Commonwealth